



# Luciferase Reporter Vector

User Manual  
Cat#: LR00XX

**Panomics, Inc.**

Luciferase Reporter Vector User Manual

**Copyright**

© Copyright 2008, Panomics, Inc. All rights reserved.

**Trademarks**

FurGene is a trademark of Hoffman Roche, OptiMem is trademark of Invitrogen. All trademarks are of their respective owners.

**Citing Luciferase Reporter Vectors in Publications**

When describing a procedure for publication using this product, we would appreciate it if you would refer to it as the Luciferase Reporter Vector from Panomics

If a paper cites the Luciferase Reporter Vector and is published in a research journal, the lead author(s) may receive a travel stipend for use at a technology conference or tradeshow by sending a copy of the paper to our technical support group at [techsupport@panomics.com](mailto:techsupport@panomics.com) or via fax at (510) 818-2610.

**Disclaimer**

Panomics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice.

Although this manual has been prepared with every precaution to ensure accuracy, Panomics, Inc. assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.

## Introduction

**Overview** Eukaryotic gene expression is regulated by a wide variety of developmental and environmental stimuli. First, an extracellular signaling molecule binds to a specific receptor. The signal is then transmitted through a series of molecular cascades, which activate or deactivate specific transcription factors (TFs) that regulate gene expression. The expression of any given gene is controlled by multiple transcription factors, which in turn are modulated by multiple signal transduction pathways. These pathways are also interconnected by molecular “cross-talk” (1–2).

In order to fully understand the mechanisms underlying gene expression, we must tease apart these multiple layers of regulation. An important step toward this goal is studying the activation of transcription factors. A popular method for this is the use of transcription reporter vectors, which contain a *cis*-acting DNA response element upstream of a reporter gene, such as luciferase. By transfecting such a vector into cells, you can examine the effects of various stimuli on a given signaling pathway *in vivo*.

### Luciferase Reporter Vectors

Panomics' Luciferase Reporter Vectors are designed to monitor transcription factor binding activity *in vivo* through the use of the standard luciferase assay. Each Luciferase Reporter Vector contains multiple repeats of a specific transcription factor binding element. Binding at this recognition site by the corresponding transcription factor results in the expression of luciferase, which initiates a powerful bioluminescent reaction when exposed to its substrate, luciferin. Light emitted from the chemical reaction is directly proportional to the amount of enzyme and thus the binding activity of the targeted transcription factor.

Currently, Luciferase Reporter Vectors are available for more than 50 transcription factors, all of which are also covered by our Protein/DNA Arrays. See [www.panomics.com](http://www.panomics.com) for an up-to-date list of all of the Luciferase Reporter Vectors

### Principle of Luciferase Reporter Vectors


The Luciferase Reporter Vectors have been specially constructed to report the binding activity of an individual TF. Multiple copies of the *cis*-acting enhancer element have been inserted into each vector upstream of a minimal TA promoter, the TATA box from the Herpes simplex virus thymidine kinase promoter. This promoter sequence drives expression of the luciferase gene (*luc*). The backbone of the vector contains an antibiotic resistance gene for cloning purposes, an origin of replication, and an *f1* origin for single-stranded DNA production.

To assess *in vivo* TF binding activity, the Luciferase Vector is first transfected into cells. If desired, an antibiotic resistant vector can be cotransfected to establish a stable cell line. After a set amount of time, the cells are lysed and subjected to the standard luciferase assay. Luminescence is detected and measured by a luminometer or scintillation counter. The resulting data can be used to quantify TF activity.

### A Valuable Tool for Monitoring Transactivation

Reporter Vectors facilitate the study of signal transduction. With the aid of the reporter vectors, *in vivo* transcription factor activity in cell lines of various origins or those treated with a stimulus of interest can easily be compared. Cotransfect a vector expressing a gene of interest along with a Luciferase Reporter Vector to observe the effects of the gene of interest on signaling pathway activity.

Reporter Vectors are also perfectly suited to be used in conjunction with Panomics' own Arrays as well. Like our popular EMSA gel-shift kits, Luciferase Reporter Vectors can be used to confirm and validate binding data from the Protein/DNA Arrays. The most popular transcription factors from the three version of the Protein/DNA Arrays can be monitored by our reporter vectors.



*The cis-acting enhancer element sequence resides between the NheI and BglII restriction sites, upstream of the TATA box promoter, which drives expression of the firefly luciferase reporter gene upon transcription factor binding. Unique restriction sites are listed on the map. The length of the enhancer element differs for each Luciferase Reporter Vector, therefore the positions indicated on the vector map should be adjusted accordingly. All of the Luciferase Vectors contain the Panomics signature sequence. These vectors are intended for research use only and should not be used commercially.*

## Materials

**Materials Provided** Upon receipt, immediately store the vectors at  $-20^{\circ}\text{C}$ .

Reporter Vector*	20 $\mu\text{L}$ of 500 ng/ $\mu\text{L}$ ; 10 $\mu\text{g}$
Control Reporter Vector	20 $\mu\text{L}$ of 500 ng/ $\mu\text{L}$ ; 10 $\mu\text{g}$

\* Sequence information is available on our website.

**Additional Items Required** Items Required but not provided

Item	Source
Transfection reagent	FuGENE™ 6 Transfection Reagent (Roche, Cat. # 1 815 019)
Reduced serum culture media	Opti-MEM™ Reduced Serum Medium (Invitrogen, Cat. # 31985-062)
Standard cell culture supplies	Major Laboratory Supplier
Luciferase reagent	Luciferase Assay System (Promega, cat.# E150C)
Luminometer	Panomics Luminometer or equivalent

**Transfection** The procedure is modified for using FuGENE 6 (Roche) in a 12-well culture plate. We do not recommend the use of other transfection methods, as this leads to reduced luciferase activity from the reporter.

*Guidelines for transfecting Luciferase Reporter Vector with FuGENE 6.*

Culture Vessel	Volume of plating medium	Vector DNA (mg) and dilution volume	FuGENE 6 Reagent
96-well	100 $\mu\text{L}$	0.05 $\mu\text{g}$ in 5 $\mu\text{L}$	0.1 - 0.3 $\mu\text{L}$
24-well	500 $\mu\text{L}$	0.2 $\mu\text{g}$ in 20 $\mu\text{L}$	0.6 - 1.8 $\mu\text{L}$
12-well	1 mL	1.0 $\mu\text{g}$ in 50 $\mu\text{L}$	1.2 - 3.6 $\mu\text{L}$
35-mm	2 mL	2.0 $\mu\text{g}$ in 100 $\mu\text{L}$	3.0 - 9.0 $\mu\text{L}$
6-well	2 mL	2.0 $\mu\text{g}$ in 100 $\mu\text{L}$	3.0 - 9.0 $\mu\text{L}$
60-mm	5 mL	5 $\mu\text{g}$ in 200 $\mu\text{L}$	6.0 -20.0 $\mu\text{L}$

1.1 The day before transfection, plate  $1-3 \times 10^5$  cells in 1 mL of their normal growth medium containing serum **without antibiotics**. Always plate cells in duplicate. This amount of cells should yield 50–80% confluence on the day of transfection.

**Note:** Lower confluence is required to allow enough surface area for growth during the experiment period.

1.2 For each well of cells to be transfected, dilute 0.5  $\mu\text{g}$  (1  $\mu\text{L}$ ) of the Luciferase

Reporter Vector or the Luciferase Control Vector with 50  $\mu$ L of Opti-MEM I Reduced Serum Medium or serum-free culture media. **Note:** Each vector should be transfected in duplicate; one group of cells is for treatment and the other group is for control.

- 1.3 For each well, add 3  $\mu$ L FuGENE 6 Reagent directly into 50  $\mu$ L of Opti-MEM I Reduced Serum Medium or serum-free culture media, mix, and incubate for 5 min at room temperature, but no longer. This dilution can be prepared in bulk for multiple wells.
- 1.4 Combine the diluted Luciferase Vector with the diluted FuGENE 6 Reagent and mix gently. Incubate for 15-45 min at room temperature to allow DNA/transfection reagent complexes to form. **Note:** Do not allow undiluted FuGENE 6 Reagent to come into contact with plastic surfaces other than pipette tips. Once the FuGENE 6 Reagent is diluted, combine it with the diluted DNA (from Step 1.2) within 45 min.
- 1.5 Add 100 mL of the DNA/transfection reagent complex directly to the complete growth media on cells and mix gently by rocking the plate back and forth.
- 1.6 Incubate the cells at 37°C in a CO<sub>2</sub> incubator, overnight.  
**Note:** It is possible to allow the transfection reaction to proceed for only 6-8 hrs. However, we recommend allowing 12 - 16 hr for gene expression, prior to treatment of cells.
- 1.7 After transfection, the media should be removed and for an untreated sample, replaced with new media; for a treated sample, replace old media with new media that has been treated with the stimulus. When removing and replacing media, use caution to ensure that the cells remain attached to the bottom of the plate.
- 1.8 Keep the transfected cells at 37°C in a CO<sub>2</sub> incubator for the remainder of the incubation period. Inducer concentrations and duration of stimulation should be optimized for cell type and inducer reagent used.

### Collecting cells

- 2.1 After incubation, aspirate to completely remove the media from the culture plates. Be careful not to disturb the cells in the process.
- 2.2 Lyse the attached cells by adding lysis buffer (Promega, Luciferase Assay System) to each well. Use approximately 50  $\mu$ L per well for a 12-well plate; 100  $\mu$ L per well for a 6-well plate.
- 2.3 To detach cells from the plate, immediately pipet the mixture up and down.
- 2.4 Transfer the cell lysate/buffer solution to a clean 1.5 mL microcentrifuge tube. Keep on ice or store at -20°C.

**Assay for luciferase activity following the instructions given by the supplier (Promega, Luciferase Assay System, cat.# E1500).**

- References**
1. Gottlicher, M., *et al.* (1998) *J Mol Med* **76**:480–489.
  2. Frønsdal, K., *et al.* (1998) *J Biol Chem* **273**: 31853–31859.
  3. Blackshear, P., *et al.* (1987) *J Biol Chem* **262**: 7774–7781.
  4. Petricoin, E., *et al.* (1992) *Mol Cell Biol* **12**: 4486 – 4495
  5. Stumpo, D., *et al.* (1988) *J Biol Chem* **263**: 1611–1614.
  6. Muller, S., *et al.* (1995) *J Invest Dermatol* **104**: 970–975.
  7. Orzechowski, H., *et al.* (2001) *Mol Pharmacol* **60**: 1332–1342.
  8. You, L and S. Jakowlew (1997) *Am J Respir Cell Mol Biol* **17**: 617–624.
  9. Suzan, M., *et al.* (1991) *J Immunol* **146**: 377 – 383.
  10. Eyries, M., *et al.* ( 2002) *Circ Res* **91**: 899–906.

## Contacting Panomics

For ordering information or technical support, contact the appropriate resource provided below according to your geographical location.

Location	U.S. Corporate Headquarters
Address	Panomics, Inc. 6519 Dumbarton Circle Fremont CA 94555 USA
Telephone	1.510.818.2600
FAX	1.510.818.2610
Email	info@panomics.com
Technical Support	<a href="mailto:techsupport@panomics.com">techsupport@panomics.com</a> or <a href="tel:1.877.726.6642">1.877.726.6642 option 3</a>
Ordering Information	orders@panomics.com

Location	European Headquarters
Address	Panomics Srl Via Sardegna 1 20060 Vignate Milano Italy
Telephone	+39.02.95.360.250
FAX	+39.02.95.360.992
Email	info_europe@panomics.com
Technical Support	<a href="mailto:techsupport_europe@panomics.com">techsupport_europe@panomics.com</a>
Ordering Information	order_europe@panomics.com

Location	Asia Pacific Headquarters
Address	Panomics, Inc. 16F Gemdale Plaza Tower A, No. 91 Jianguo Road Beijing 100022 P.R. China
Telephone	+86.10.59208157
FAX	+86.10.59208111
Email	info_asia@panomics.com
Technical Support	<a href="mailto:techsupport_asia@panomics.com">techsupport_asia@panomics.com</a>
Ordering Information	order_asia@panomics.com

For an updated list of FAQs and product support literature, visit our website at [www.panomics.com](http://www.panomics.com).

---