



Telephone: 1-866-757-2414 (Toll Free) 604-247-2416 (Local)
Address: #8-13520 Crestwood Place Richmond, BC, Canada V6V 2G2

iLenti™-GFP Scramble siRNA

Catalog Number LV015-G

Original Vector	Cat#LV016 – Bbs I Linearized iLenti™ GFP siRNA Vector
Cloning Site	BbsI (cut site on the vector is eliminated after cloning)
Promoter	hH1 and hU6
Insert	5' – GGGTGAAGTACGTCAGAA – 3' [Negative Control]
Insert Reference	Yu EZ, Li YY, Liu XH, Kagan E, McCarron RM. 2004. Lab Invest. 84(5):553-61.
Tag	None
Reporter	GFP
Selection Marker	Bacterial System: Kanamycin Mammalian System: Puromycin
Amount	500ng
Appearance	Liquid
Storage Temperature	-20°C or below (shipped at ambient temperature)
Shelf Life	1 year (when at -20°C or below in a non-frost free freezer)
QC Method	Restriction Digestion and Sequencing
Sequencing Primer	5' – CAACCCGCTCCAAGGAATCG – 3'

