









NeuroPORTER Transfection Data





C6 Glioma Cells
transfected with 30 μ l
NeuroPORTER and 3 μ g of
gWiz/GFP vector.



Rat Primary Cortical Neurons
transfected with 10 μ l
NeuroPORTER and 2 μ g of
gWiz/GFP vector.



Human NT2 Cells
transfected with 5 μ l
NeuroPORTER and 2 μ g of
gWiz/GFP vector.



Differentiated hNT Neurons
transfected with 20 μ l
NeuroPORTER and 2 μ g of
gWiz/GFP vector.

Gene Therapy Systems, Inc.
10190 Telesis Court, San Diego, Ca 92121
Voice: (858) 457-1919 • 888-428-0558
Fax Lines: (858) 623-9494 • (858) 587-1499
website: <http://www.genetherapysystems.com>