



The New Standard in Exosome Purification

Introducing the Most Comprehensive Range of Exosome Purification Kits Available Anywhere


- ◆ **Plasma/Serum:** purify exosomes from 50 μ L to 10 mL inputs
- ◆ **Urine:** purify exosomes from 250 μ L to 30 mL inputs
- ◆ **Rapid and Simple:** purify and enrich exosomes in under 45 minutes
- ◆ **Quality:** as good as ultracentrifugation
- ◆ **Convenient:** no filtration, no tedious centrifugation, precipitation reagents or special syringes required
- ◆ **Applications:** Exosomes suitable for Nanosight® analysis, EM and Western Blot
Purify exosome-specific miRNA for NGS, qRT-PCR, microarrays, gene expression analysis, biomarker discovery and more...




Highest Concentration of Exosomes



Rapid Isolation of Pure, Intact Exosomes



Typical Size Range of Exosomes & Extracellular Vesicles



Excellent Linearity for Exosomal RNA Recovery

The New Standard in Exosome Purification

Ordering Information

Plasma Exosome Kits	Cat.#	Size (Preps)	Input Range
Plasma/Serum Exosome Purification Mini Kit	57400	50	50 μ L to 1 mL
Plasma/Serum Exosome Purification Midi Kit	57500	25	1 mL to 4 mL
Plasma/Serum Exosome Purification Maxi Kit	57600	15	4 mL to 10 mL

Urine Exosome Kits	Cat.#	Size (Preps)	Input Range
Urine Exosome Purification Mini Kit	57700	50	250 μ L to 1 mL
Urine Exosome Purification Midi Kit	57800	25	2 mL to 10 mL
Urine Exosome Purification Maxi Kit	57900	15	11 mL to 30 mL

Exosome Kit	Cat.#	Size (Preps)
Exosomal RNA Isolation Kit	58000	50

Exosome Isolation and RNA Purification Kit	Cat.#	Size (Preps)	Input Range
Plasma/Serum Exosome Purification and RNA Isolation Mini Kit	58300	50	50 μ L to 1 mL
Plasma/Serum Exosome Purification and RNA Isolation Midi Kit	58500	25	1 mL to 4 mL
Plasma/Serum Exosome Purification and RNA Isolation Maxi Kit	58600	15	4 mL to 10 mL
Urine Exosome Purification and RNA Isolation Mini Kit	58400	50	250 μ L to 1 mL
Urine Exosome Purification and RNA Isolation Midi Kit	58700	25	2 mL to 10 mL
Urine Exosome Purification and RNA Isolation Maxi Kit	58800	15	11 mL to 30 mL