


# NEW: Preserved Blood RNA Purification Kits

For rapid and simple isolation of total RNA from blood preserved in Tempus™ or PAXgene™ Blood RNA Tubes


## Features:

- ◆ Isolate all sizes of RNA including microRNA without phenol
- ◆ Choose optimized protocols for either Tempus™ or PAXgene™ blood tubes
- ◆ Purified RNA is of excellent quality and ready for any downstream application
- ◆ Rapid and convenient spin column protocols




### Isolate a Diversity of RNA Species

Total RNA was isolated from 3 mL hamster blood samples collected into Tempus™ RNA Blood Tubes using Norgen's Preserved Blood RNA Purification Kit I. The purified total RNA was then resolved on an Agilent® BioAnalyzer. A: Electropherogram, B: Gel Banding Pattern.




### RNA Ready for Any Downstream Application

Total RNA was isolated from human blood collected in a PAXgene™ Blood RNA Tube followed by RNA isolation using Norgen's Preserved Blood RNA Purification Kit II. The purified RNA was then used as the template in a RT-qPCR for detecting the S15 gene and for detecting miR-21. Norgen's kit isolated complete range of RNA including high yields of mRNA and microRNA.



Norgen Competitor



Small RNA Species

**Isolate High Quality RNA with Complete Size Range.** Unlike most competitor's kits, Norgen's Preserved Blood RNA Purification Kits allow for the isolation of all sizes of RNA, from the very large RNA down to the microRNA without the use of phenol.

*Tempus™ is a registered trademark of Life Technologies  
PAXgene™ is a registered trademark of PreAnalytix GmbH*

**NORGEN**  
BIOTEK  CORP.

## Ordering Information

Kit	Catalogue Number
Preserved Blood RNA Purification Kit I (for use with Tempus™ Blood RNA Tubes)	43400
Preserved Blood RNA Purification Kit II (for use with PAXgene™ Blood RNA Tubes)	43500

### PLEASE CONTACT US FOR MORE INFORMATION:

North American Toll Free: 1-866-667-4362  
Phone: 905-227-8848  
Fax: 905-227-1061

[www.norgenbiotech.com](http://www.norgenbiotech.com)  
[info@norgenbiotech.com](mailto:info@norgenbiotech.com)

3430 Schmon Parkway, Thorold, ON, Canada L2V 4Y6  
ISO 13485:2003, ISO 9001:2008 Registered

Ver1.0