

Product name **3-D Life Thioglycerol**

Description Thioglycerol (CAS-No. 96-27-5) can be covalently coupled with its thiol group to SH-reactive groups of Dextran or PVA polymers to modify 3-D Life Hydrogels.

Catalog number T10-3

Quantity Allows modification of 3 ml 3-D Life Hydrogel with 0.4 mMol/L thioglycerol.

Applications 3-D Life Thioglycerol can be used as a negative control in place of peptides to modify hydrogels.

Components

	Materials provided ¹	Quantity	Concentration reactive groups	Storage Temperature
	Thioglycerol	180 µl	20 mmol/l	-20°C or -80°C

¹ All materials are filter-sterilized.

NOTE INTENDED FOR RESEARCH USE ONLY. NOT FOR USE IN HUMAN THERAPEUTIC OR DIAGNOSTIC APPLICATIONS.